

Product datasheet for **HP231699**

RBIS Human qPCR Primer Pair (NM_001099671)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	401466
Forward Sequence:	AAAGCAAACCCAGTTACCACTAATC
Reverse Sequence:	GGTTCAAGTGAAATGCTTTTTGCG
ACCN:	BC032347 , NM_001099671 , NM_001099671.1 , NM_001099671.2 , BC032347.1 , BC055092 , BC062221 , BC146983 , BC146988 , BC157825 , BC171873 , BG926519 , BM785759
UniProt ID:	Q8N0T1
Synonyms:	C8orf59
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)