

## Product datasheet for **HP231702**

### DIP2 homolog B (DIP2B) Human qPCR Primer Pair (NM\_173602)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	57609
Forward Sequence:	GTTGCGGAACAAAGACCTGATGC
Reverse Sequence:	ATTGGCTGGCACCAGAGCAAGA
ACCN:	<a href="#">NM_173602</a> , <a href="#">NM_173602.2</a> , <a href="#">BC007671</a> , <a href="#">BC030156</a> , <a href="#">BC075027</a> , <a href="#">BU101773</a> , <a href="#">BU102128</a> , <a href="#">BU429622</a> , <a href="#">BX537905</a> , <a href="#">NM_173602.3</a>
UniProt ID:	<a href="#">Q9P265</a>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)