

Product datasheet for **HP231717**

TAX1BP1 Human qPCR Primer Pair (NM_006024)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8887
Forward Sequence:	GTTAGCTGATGCAGTGGCAGAAC
Reverse Sequence:	CAGCCATCTGAAGACGGAGTTTC
ACCN:	NM_006024 , NM_006024.1 , NM_006024.2 , NM_006024.3 , NM_006024.4 , NM_006024.5 , NM_006024.6 , BC050358 , BC050358.2 , BC024600 , BP280714 , BU679319
UniProt ID:	Q86VP1
Synonyms:	CALCOCO3; T6BP; TXBP151
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)