

Product datasheet for **HP231750**

MTX2 Human qPCR Primer Pair (NM_006554)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10651
Forward Sequence:	GTATGGATCTCCTTACCCTTGCC
Reverse Sequence:	CTTGACAGCACTGGTCTACATCC
ACCN:	NM_006554 , NM_006554.1 , NM_006554.2 , NM_006554.3 , NM_006554.4 , BC017271 , BC017271.2 , BC067831 , BC067831.1 , BC088359 , BG699784 , BU165147 , NM_006554.5
UniProt ID:	O75431
Synonyms:	MDPS; metaxin-2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)