

## Product datasheet for **HP231780**

### ZSCAN23 Human qPCR Primer Pair (NM\_001012455)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	222696
Forward Sequence:	GAGAAGGCAACTCCAGGAGCAG
Reverse Sequence:	CCACATTCCAAGTCCCAGCCTT
ACCN:	<a href="#">NM_001012455</a> , <a href="#">NM_001012455.1</a> , <a href="#">BC101557</a> , <a href="#">BC101559</a> , <a href="#">BC143560</a> , <a href="#">BX283350</a> , <a href="#">NM_001012455.2</a>
UniProt ID:	<a href="#">Q3MJ62</a>
Synonyms:	dj29K1.3; dj29K1.3.1; ZNF390; ZNF453
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)