

Product datasheet for **HP231789**

MTGR1 (CBFA2T2) Human qPCR Primer Pair (NM_005093)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9139
Forward Sequence:	GGAGCAAACCATAGCGGATGTC
Reverse Sequence:	GGCAGAAAGAGCCACAGTATCG
ACCN:	NM_005093 , NM_005093.1 , NM_005093.2 , NM_005093.3 , BC040344 , BC040344.2 , BC015066 , BC016298
UniProt ID:	O43439
Synonyms:	EHT; MTGR1; p85; ZMYND3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)