

Product datasheet for **HP231790**

ZGPAT Human qPCR Primer Pair (NM_032527)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	84619
Forward Sequence:	CTCCTCACCAAGATGGGCTATG
Reverse Sequence:	TTGCCAACCTGGTCTGCTTCT
ACCN:	NM_032527 , NM_032527.1 , NM_032527.2 , NM_032527.3 , NM_032527.4 , BC032612 , BC032612.1 , BC019338
UniProt ID:	Q8N5A5
Synonyms:	GPATC6; GPATCH6; KIAA1847; ZC3H9; ZC3HDC9; ZIP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)