

Product datasheet for **HP231835**

Leptin Receptor (LEPR) Human qPCR Primer Pair (NM_002303)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3953
Forward Sequence:	GCAGTCTATGCTGTTTCAGGTGC
Reverse Sequence:	CCAAAATTCAGGTCCTCTCATAGG
ACCN:	NM_002303 , NM_002303.1 , NM_002303.2 , NM_002303.3 , NM_002303.4 , NM_002303.5 , BC031086 , BC035120 , BC131779 , BC160009 , BX475249 , NM_002303.6
UniProt ID:	P48357
Synonyms:	CD295; LEP-R; LEPRD; OB-R; OBR
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)