

Product datasheet for **HP231839**

Aprataxin (APTX) Human qPCR Primer Pair (NM_175073)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	54840
Forward Sequence:	GGAAGCAGTTGTGATTGGGCGT
Reverse Sequence:	ACGACTGAGTCAATGCTGGTGG
ACCN:	NM_175073 , NM_175073.1 , NM_175073.2 , BC032650 , BC001628 , BC068107 , BC104881 , BQ957536 , BX538161
UniProt ID:	Q7Z2E3
Synonyms:	AOA; AOA1; AXA1; EAOH; EOAHA; FHA-HIT
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)