

## Product datasheet for **HP231853**

### SHC (SHC1) Human qPCR Primer Pair (NM\_183001)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6464
Forward Sequence:	ACAGCCGAGTATGTCGCCTATG
Reverse Sequence:	CAATGGTGCTGATGACATCCTGG
ACCN:	<a href="#">NM_183001</a> , <a href="#">NM_183001.1</a> , <a href="#">NM_183001.2</a> , <a href="#">NM_183001.3</a> , <a href="#">NM_183001.4</a> , <a href="#">BC014158</a> , <a href="#">BC014158.2</a> , <a href="#">BC011813</a> , <a href="#">BC017283</a> , <a href="#">BC033925</a> , <a href="#">BX647136</a> , <a href="#">BX647149</a> , <a href="#">NM_183001.5</a>
UniProt ID:	<a href="#">P29353</a>
Synonyms:	SHC; SHCA
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)