

Product datasheet for **HP231929**

CYLD Human qPCR Primer Pair (NM_015247)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1540
Forward Sequence:	GGTAATCCGTTGGATCGGTCAG
Reverse Sequence:	AGTGCCTCTGAAGTTCCATCC
ACCN:	<u>NM_015247</u> , <u>NM_015247.1</u> , <u>NM_015247.2</u> , <u>BC012342</u> , <u>BC012342.1</u> , <u>BI458962</u> , <u>BX642061</u> , <u>NM_015247.3</u>
UniProt ID:	<u>Q9NQC7</u>
Synonyms:	BRSS; CDMT; CYLD1; CYLDI; EAC; FTDALS8; MFT; MFT1; SBS; TEM; USPL2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)