

Product datasheet for **HP231976**

SNX29 Human qPCR Primer Pair (NM_001080530)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	92017
Forward Sequence:	GAAGTCCTCTGAGCAGCCTGTT
Reverse Sequence:	CGCAGTGATCTGTTCTCCTCT
ACCN:	NM_001080530 , NM_001080530.1 , BC013346 , BC029857 , BC156063 , BC172456 , BM692924 , BQ185857 , BX648908
UniProt ID:	Q8TEQ0
Synonyms:	A-388D4.1; FLJ00143
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)