

Product datasheet for **HP231984**

CSNK1G1 Human qPCR Primer Pair (NM_022048)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	53944
Forward Sequence:	GGAGGACTTGTTTGACCTCTGTG
Reverse Sequence:	GACCAATCAGGAAGTTCTCTGGC
ACCN:	<u>NM_022048</u> , <u>NM_022048.1</u> , <u>NM_022048.3</u> , <u>NM_022048.4</u> , <u>BC130613</u> , <u>BC017236</u> , <u>BC062995</u> , <u>BC077731</u> , <u>BC136504</u> , <u>BE644696</u> , <u>NM_022048.5</u>
UniProt ID:	<u>Q9HCP0</u>
Synonyms:	CK1gamma1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)