

Product datasheet for **HP232082**

Neuronal calcium binding protein (NECAB3) Human qPCR Primer Pair (NM_031232)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	63941
Forward Sequence:	GCATTGGAATCGCTGAACCGTG
Reverse Sequence:	CGCGTCACAAACTGGTCCACTT
ACCN:	NM_031232 , NM_031232.1 , NM_031232.2 , NM_031232.3 , BC047673 , BC023270 , BQ212085 , BQ431858 , NM_031232.4
UniProt ID:	Q96P71
Synonyms:	APBA2BP; dj63M2.4; dj63M2.5; EFCBP3; NIP1; STIP3; SYTIP2; XB51
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)