

Product datasheet for **HP232136**

POM121 Human qPCR Primer Pair (NM_172020)

Product data:

| | |
|-------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 9883 |
| Forward Sequence: | ATCACTGCCGAGGACCTAGACT |
| Reverse Sequence: | GGTCTCAGTGACAGATTCGAG |
| ACCN: | NM_172020 , NM_172020.1 , NM_172020.2 , NM_172020.3 , NM_172020.4 , BC008794 , BC008794.2 , BC001518 , BE408761 , BG831530 , BI544820 , BM994161 |
| UniProt ID: | Q96HA1 |
| Synonyms: | P145; POM121A |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions) |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)