

Product datasheet for **HP232149**

CPSF30 (CPSF4) Human qPCR Primer Pair (NM_006693)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10898
Forward Sequence:	ACCAAGATGCCCCGAGTGCTACT
Reverse Sequence:	CCACGGTCATACCAAGGACAGT
ACCN:	<u>BC050738</u> , <u>NM_006693</u> , <u>NM_006693.1</u> , <u>NM_006693.2</u> , <u>NM_006693.3</u> , <u>BC003101</u> , <u>BC003101.1</u> , <u>BE388179</u> , <u>BM757197</u> , <u>BM787372</u> , <u>BQ052009</u> , <u>BQ641871</u> , <u>BX358595</u>
UniProt ID:	<u>O95639</u>
Synonyms:	CPSF30; NAR; NEB-1; NEB1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)