

## Product datasheet for **HP232187**

### COMMD9 Human qPCR Primer Pair (NM\_014186)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	29099
Forward Sequence:	CATGTTCCAGCTTGTCTGTGACC
Reverse Sequence:	AAGAGAGCCAGAATTGCCTCGG
ACCN:	<a href="#">NM_014186</a> , <a href="#">NM_014186.1</a> , <a href="#">NM_014186.2</a> , <a href="#">NM_014186.3</a> , <a href="#">BC010892</a> , <a href="#">BC010892.1</a> , <a href="#">BC092458</a> , <a href="#">BC108685</a> , <a href="#">BG702189</a> , <a href="#">BI668465</a> , <a href="#">BM172201</a> , <a href="#">NM_014186.4</a>
UniProt ID:	<a href="#">Q9P000</a>
Synonyms:	C11orf55; HSPC166; LINC00610
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)