

Product datasheet for **HP232228**

BPTF Human qPCR Primer Pair (NM_182641)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2186
Forward Sequence:	GGAGAGATGTTGGTCCTTATGCC
Reverse Sequence:	CTTCCTCTGAGGTGTAGCGT
ACCN:	<u>NM_182641</u> , <u>NM_182641.1</u> , <u>NM_182641.2</u> , <u>NM_182641.3</u> , <u>BC037422</u> , <u>BC039881</u> , <u>BC050566</u> , <u>BC067234</u> , <u>BE464864</u> , <u>BQ716216</u> , <u>BU169920</u> , <u>BU176205</u> , <u>NM_182641.4</u>
UniProt ID:	<u>Q12830</u>
Synonyms:	FAC1; FALZ; NEDDFL; NURF301
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)