

## Product datasheet for **HP232234**

### GPR18 Human qPCR Primer Pair (NM\_005292)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2841
Forward Sequence:	CTTCTGCCAGATTCTTGGAGCTC
Reverse Sequence:	GTTCTTTGGCGTACTTCGGCTG
ACCN:	<a href="#">NM_005292</a> , <a href="#">NM_005292.1</a> , <a href="#">NM_005292.2</a> , <a href="#">NM_005292.3</a> , <a href="#">BC050646</a> , <a href="#">BC050646.1</a> , <a href="#">BC008569</a> , <a href="#">BC066927</a> , <a href="#">BI770379</a> , <a href="#">BT009907</a>
UniProt ID:	<a href="#">Q14330</a>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)