

Product datasheet for **HP232270**

TCP1 eta (CCT7) Human qPCR Primer Pair (NM_006429)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10574
Forward Sequence:	TGCCAGGTGTTTGAAGAGACCC
Reverse Sequence:	GTCTCCTCCATAAACTGCTCGG
ACCN:	NM_006429 , NM_006429.1 , NM_006429.2 , NM_006429.3 , BC088351 , BC088351.1 , BC019296 , BM558355 , BM918058 , NM_006429.4
UniProt ID:	Q99832
Synonyms:	CCTETA; CCTH; NIP7-1; TCP1ETA
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)