

Product datasheet for **HP232272**

CrkRS (CDK12) Human qPCR Primer Pair (NM_016507)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	51755
Forward Sequence:	CGAACTCAGCAAAATGGCTCCTC
Reverse Sequence:	TTGGATGGAGGTGGCTCTTCGA
ACCN:	NM_016507 , NM_016507.1 , NM_016507.2 , NM_016507.3 , BC140854 , BC150265 , BE003756 , BM557102 , NM_016507.4
UniProt ID:	Q9NYV4
Synonyms:	CRK7; CRKR; CRKRS
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)