

Product datasheet for **HP232281**

STEAP3 Human qPCR Primer Pair (NM_018234)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55240
Forward Sequence:	TGCAAACCTCGCTCAACTGGAGG
Reverse Sequence:	AGGCAGGTAGAACTTGTAGCGG
ACCN:	<u>NM_018234</u> , <u>NM_018234.1</u> , <u>NM_018234.2</u> , <u>BC095421</u> , <u>BC095421.1</u> , <u>BC042150</u> , <u>BM146930</u> , <u>BX538047</u> , <u>NM_018234.3</u>
UniProt ID:	<u>Q658P3</u>
Synonyms:	AHMIO2; dudlin-2; dudulin-2; pHyde; STMP3; TSAP6
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)