

Product datasheet for **HP232315**

ENPP2 Human qPCR Primer Pair (NM_006209)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5168
Forward Sequence:	TATGCTGCGGAAACTCGTCAGG
Reverse Sequence:	GACGTTGACACACCGATGCAGT
ACCN:	NM_006209 , NM_006209.1 , NM_006209.2 , NM_006209.3 , NM_006209.4 , BC034961 , BM694754 , BQ365938 , BU731264 , NM_006209.5
UniProt ID:	Q13822
Synonyms:	ATX; ATX-X; AUTOTAXIN; LysoPLD; NPP2; PD-IALPHA; PDNP2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)