

Product datasheet for **HP232379**

NFAT4 (NFATC3) Human qPCR Primer Pair (NM_173165)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4775
Forward Sequence:	AGACAGTCGCTACTGCAAGCCA
Reverse Sequence:	GCGGAGTTTCAAATACCTGCAC
ACCN:	<u>NM_173165</u> , <u>NM_173165.1</u> , <u>NM_173165.2</u> , <u>BC001050</u> , <u>BC001050.2</u> , <u>BE350846</u> , <u>BU849493</u> , <u>NM_173165.3</u>
UniProt ID:	<u>Q12968</u>
Synonyms:	NF-AT4c; NFAT4; NFATX
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)