

Product datasheet for **HP232392**

CPEB1 Human qPCR Primer Pair (NM_030594)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	64506
Forward Sequence:	TGGCAGCCATCTTGAACGACCT
Reverse Sequence:	GCTGACTGCTTTCAGGTAAGTCC
ACCN:	<u>BC035348</u> , <u>NM_030594</u> , <u>NM_030594.1</u> , <u>NM_030594.2</u> , <u>NM_030594.3</u> , <u>NM_030594.4</u> , <u>BC035348.1</u> , <u>BC050629</u> , <u>BC050629.1</u> , <u>BM715948</u> , <u>BX649159</u> , <u>NM_030594.5</u>
UniProt ID:	<u>Q9BZB8</u>
Synonyms:	CPE-BP1; CPEB; CPEB-1; h-CPEB; hCPEB-1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)