

Product datasheet for **HP232441**

C3orf55 (PQLC2L) Human qPCR Primer Pair (NM_001130002)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	152078
Forward Sequence:	GAGCACTTGACCTGCCTTAGAAG
Reverse Sequence:	GGCAGTTGGTTAGTCAGGTAGC
ACCN:	<u>NM_001130002</u> , <u>NM_001130002.1</u> , <u>NM_001130002.2</u> , <u>BC010062</u> , <u>BC029129</u> , <u>BC126187</u> , <u>BC126189</u> , <u>BC144066</u> , <u>BE905784</u> , <u>BG209594</u> , <u>BG217991</u> , <u>BX113701</u>
UniProt ID:	<u>A1A4F0</u>
Synonyms:	C3orf55; PQLC2L; SLC66A2L
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)