

Product datasheet for **HP232450**

TRPM3 Human qPCR Primer Pair (NM_206946)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	80036
Forward Sequence:	GGAAGTGTTTGC GGACCAGATAG
Reverse Sequence:	GCACGATCCAAGCTCCTGTCTT
ACCN:	<u>NM_206946</u> , <u>NM_206946.1</u> , <u>NM_206946.2</u> , <u>NM_206946.3</u> , <u>BC022454</u> , <u>BC067733</u> , <u>BC094699</u> , <u>BC121821</u> , <u>BC134414</u> , <u>BC142972</u> , <u>BC172350</u> , <u>BM668925</u> , <u>NM_206946.5</u>
UniProt ID:	<u>Q9HCF6</u>
Synonyms:	GON-2; LTRPC3; MLSN2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)