

Product datasheet for **HP232458**

CTIP1 (BCL11A) Human qPCR Primer Pair (NM_022893)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	53335
Forward Sequence:	CTCGTTCTGCACATGGAGCTCT
Reverse Sequence:	GCAAGAGAAACCATGCACTGGTG
ACCN:	<u>NM_022893</u> , <u>NM_022893.1</u> , <u>NM_022893.2</u> , <u>NM_022893.3</u> , <u>BC021098</u> , <u>BC021098.1</u> , <u>BG393787</u>
UniProt ID:	<u>Q9H165</u>
Synonyms:	BCL11A-L; BCL11a-M; BCL11A-S; BCL11A-XL; CTIP1; DILOS; EVI9; HBFQTL5; ZNF856
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)