

## Product datasheet for **HP232476**

### KCNC4 Human qPCR Primer Pair (NM\_004978)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3749
Forward Sequence:	CCATGTACTGCAAGTCTGAGGAG
Reverse Sequence:	CCTCATCAGACAGCACTGCGTT
ACCN:	<a href="#">NM_004978</a> , <a href="#">NM_004978.1</a> , <a href="#">NM_004978.2</a> , <a href="#">NM_004978.3</a> , <a href="#">NM_004978.4</a> , <a href="#">BC101769</a> , <a href="#">BC101769.1</a> , <a href="#">BC025377</a> , <a href="#">BF969982</a> , <a href="#">BM765446</a> , <a href="#">BX098184</a> , <a href="#">BX537638</a> , <a href="#">NM_004978.6</a>
UniProt ID:	<a href="#">Q03721</a>
Synonyms:	C1orf30; HKSHIIC; KSHIIC; KV3.4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)