

Product datasheet for **HP232569**

DcR3 (TNFRSF6B) Human qPCR Primer Pair (NM_032945)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8771
Forward Sequence:	CTCAATGTGCCAGGCTCTTCCT
Reverse Sequence:	GGAGATGTCCTGGAAAGCCACA
ACCN:	<u>NM_032945</u> , <u>NM_032945.1</u> , <u>NM_032945.2</u> , <u>BC017065</u> , <u>BC017065.1</u> , <u>BC034349</u> , <u>BM759033</u>
UniProt ID:	<u>O95407</u>
Synonyms:	DCR3; decoy receptor 3; DJ583P15.1.1; M68; OTTHUMP00000031583; TR6; tumor necrosis factor receptor superfamily, member 6b; tumor necrosis factor receptor superfamily, member 6b, decoy
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)