

## Product datasheet for **HP232639**

### SDHC Human qPCR Primer Pair (NM\_003001)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6391
Forward Sequence:	GGTCAAACCGTCCTCTGTCTC
Reverse Sequence:	CGACATGCCAAAAAGAGACCC
ACCN:	<a href="#">NM_003001</a> , <a href="#">NM_003001.1</a> , <a href="#">NM_003001.2</a> , <a href="#">NM_003001.3</a> , <a href="#">BC020808</a> , <a href="#">BC020808.1</a> , <a href="#">BC012735</a> , <a href="#">BC033626</a> , <a href="#">BC066329</a> , <a href="#">BG724095</a> , <a href="#">NM_003001.5</a>
UniProt ID:	<a href="#">Q99643</a>
Synonyms:	CYB560; CYBL; PGL3; QPS1; SDH3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)