

Product datasheet for **HP232653**

DGKG Human qPCR Primer Pair (NM_001346)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1608
Forward Sequence:	CTCCAGACTTCCGTGTTTTGGC
Reverse Sequence:	GGTCATTTCTGTCCAAGAGGC
ACCN:	<u>NM_001346</u> , <u>NM_001346.1</u> , <u>NM_001346.2</u> , <u>BC112363</u> , <u>BC089411</u> , <u>BM669549</u> , <u>NM_001346.3</u>
UniProt ID:	<u>P49619</u>
Synonyms:	DAGK3; DGK-GAMMA
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)