

Product datasheet for **HP232675**

LRRC8B Human qPCR Primer Pair (NM_001134476)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	23507
Forward Sequence:	CCACTTGACCTTCATTCCAGAAG
Reverse Sequence:	CTGAAACAGCCCATCTGGTAGC
ACCN:	<u>NM_001134476</u> , <u>NM_001134476.1</u> , <u>BC030607</u> , <u>BC030607.2</u> , <u>BC060782</u> , <u>BC035081</u> , <u>BC037831</u> , <u>BC053565</u> , <u>BI459217</u> , <u>NM_001134476.2</u>
UniProt ID:	<u>Q6P9F7</u>
Synonyms:	TA-LRRP; TALRRP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)