

Product datasheet for **HP232702**

EYA2 Human qPCR Primer Pair (NM_172113)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2139
Forward Sequence:	TCCTCCAGGAGGCATCTCACAA
Reverse Sequence:	TCTGTGTCTCCCTCTTTTCGCTG
ACCN:	<u>NM_172113</u> , <u>NM_172113.1</u> , <u>BC008803</u> , <u>BC008803.2</u> , <u>BC000289</u> , <u>BC013882</u> , <u>BT006682</u>
UniProt ID:	<u>O00167</u>
Synonyms:	EAB1; eyes absent 2; eyes absent homolog 2 (Drosophila); MGC10614; OTTHUMP00000031666
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)