

Product datasheet for **HP232722**

ZWINT Human qPCR Primer Pair (NM_007057)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	11130
Forward Sequence:	AGGCATCTTGGAACCTGTAGGC
Reverse Sequence:	GGAAATCCGCTACCTGAAGCTG
ACCN:	<u>NM_007057</u> , <u>NM_007057.1</u> , <u>NM_007057.2</u> , <u>NM_007057.3</u> , <u>BC110399</u> , <u>BC110399.1</u> , <u>BC000411</u> , <u>BC016357</u> , <u>BC020979</u> , <u>BG615304</u> , <u>BM724168</u>
UniProt ID:	<u>O95229</u>
Synonyms:	HZwint-1; KNTC2AP; SIP30; ZWINT1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)