

Product datasheet for **HP232728**

RAP1GDS1 Human qPCR Primer Pair (NM_021159)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5910
Forward Sequence:	CCACAGTCTTCCTGCAAAGCCA
Reverse Sequence:	AGCAGCACTTCCTGGTCTTTGC
ACCN:	NM_021159 , NM_021159.1 , NM_021159.2 , NM_021159.3 , NM_021159.4 , BC001816 , BC001816.2 , BC001851 , BC098269 , BC098334 , BC099708 , BC099845 , BT006837
UniProt ID:	P52306
Synonyms:	GDS1; SmgGDS
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)