

Product datasheet for **HP232865**

AGXT2L2 (PHYKPL) Human qPCR Primer Pair (NM_153373)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	85007
Forward Sequence:	GACACCAGCAACTGAAGAGGCT
Reverse Sequence:	GCATTGTCCAGGCTGAAGCACA
ACCN:	NM_153373 , NM_153373.1 , NM_153373.2 , NM_153373.3 , BC037567 , BC037567.1 , BC008009 , BC028002 , BC037238 , BC039012 , BC067745 , BC088800 , BC110335 , BX417194 , NM_153373.4
UniProt ID:	Q8IUZ5
Synonyms:	AGXT2L2; PHLU
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)