

## Product datasheet for **HP232966**

### TMEM91 Human qPCR Primer Pair (NM\_001098821)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	641649
Forward Sequence:	CAGTTCCTGTCACCGCCTCTTC
Reverse Sequence:	GTGAAAGACGGCTGCCTCCATC
ACCN:	<a href="#">NM_001098821</a> , <a href="#">NM_001098821.1</a> , <a href="#">BC063705</a> , <a href="#">BF684543</a> , <a href="#">BG697182</a> , <a href="#">BM675470</a> , <a href="#">BM726200</a> , <a href="#">BM931402</a> , <a href="#">BX090160</a> , <a href="#">NM_001098821.2</a>
UniProt ID:	<a href="#">Q6ZNR0</a>
Synonyms:	DSPC3; IFITMD6
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)