

## Product datasheet for **HP232990**

### Dynein heavy chain (DNAH9) Human qPCR Primer Pair (NM\_004662)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1770
Forward Sequence:	CCAGGCAACATGGACTACAATGG
Reverse Sequence:	CAGGAAGCCAATCTCTGCGTTC
ACCN:	<u><a href="#">NM_004662</a></u> , <u><a href="#">NM_004662.1</a></u> , <u><a href="#">NM_004662.2</a></u> , <u><a href="#">BC128261</a></u> , <u><a href="#">BC128421</a></u> , <u><a href="#">BC153883</a></u>
UniProt ID:	<u><a href="#">Q9NYC9</a></u>
Synonyms:	CILD40; DNAH17L; Dnahc9; DNEL1; DYH9; HL-20; HL20
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)