

## Product datasheet for **HP233023**

### NFAT1 (NFATC2) Human qPCR Primer Pair (NM\_173091)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4773
Forward Sequence:	GATAGTGGGCAACACCAAAGTCC
Reverse Sequence:	TCTCGCCTTCCGCAGCTCAAT
ACCN:	<u><a href="#">NM_173091</a></u> , <u><a href="#">NM_173091.1</a></u> , <u><a href="#">NM_173091.2</a></u> , <u><a href="#">NM_173091.3</a></u> , <u><a href="#">BC136418</a></u> , <u><a href="#">BC144074</a></u> , <u><a href="#">BE674599</a></u> , <u><a href="#">BF914609</a></u> , <u><a href="#">BM463165</a></u> , <u><a href="#">BU429325</a></u> , <u><a href="#">NM_173091.4</a></u>
UniProt ID:	<u><a href="#">Q13469</a></u>
Synonyms:	NFAT1; NFATP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)