

Product datasheet for **HP233027**

WDR1 Human qPCR Primer Pair (NM_017491)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9948
Forward Sequence:	AGAAGGACCACCTGCTCAGTGT
Reverse Sequence:	ATGCACCGTCAGACACTGGATC
ACCN:	NM_017491 , NM_017491.1 , NM_017491.2 , NM_017491.3 , BC002489 , BC002489.2 , BC000201 , BC030541 , BG724171 , BI831941 , BX648190 , NM_017491.5
UniProt ID:	O75083
Synonyms:	AIP1; HEL-S-52; NORI-1; PFITS
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)