

Product datasheet for **HP233142**

CD20 (MS4A1) Human qPCR Primer Pair (NM_152866)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	931
Forward Sequence:	CTGGTCCAAAACCACTCTTCAGG
Reverse Sequence:	GGCAATGTGGAAGAGCCCATTC
ACCN:	<u>NM_152866, NM_152866.1, NM_152866.2, BC002807, BC002807.2, NM_152866.3</u>
UniProt ID:	<u>P11836</u>
Synonyms:	B1; Bp35; CD20; CVID5; FMC7; LEU-16; MS4A2; S7
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)