

## Product datasheet for **HP233159**

### EML1 Human qPCR Primer Pair (NM\_004434)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2009
Forward Sequence:	GGAGGAACCAATCTCTGTGCTG
Reverse Sequence:	GGTGAAATCCGCAGCAAACAC
ACCN:	<a href="#">NM_004434</a> , <a href="#">NM_004434.2</a> , <a href="#">BC033043</a> , <a href="#">BC033043.1</a> , <a href="#">BC032541</a> , <a href="#">BX247979</a> , <a href="#">BX420178</a> , <a href="#">BX442035</a>
UniProt ID:	<a href="#">O00423</a>
Synonyms:	BH; ELP79; EMAP; EMAP-1; EMAPL
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)