

Product datasheet for **HP233192**

GGT7 Human qPCR Primer Pair (NM_178026)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2686
Forward Sequence:	GGACTTCAGCAATTACAGCGCC
Reverse Sequence:	CTCCAGGATGTTGAGAGCACTG
ACCN:	NM_178026 , NM_178026.1 , NM_178026.2 , BC141843 , BC033745 , BC040030 , BC101370 , BC101371 , BC101372 , BQ013461 , BX401449 , NM_178026.3
UniProt ID:	Q9UJ14
Synonyms:	D20S101; GGT4; GGTL3; GGTL5
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)