

Product datasheet for **HP233245**

RABL2B Human qPCR Primer Pair (NM_001130923)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	11158
Forward Sequence:	ACGTGAAGATCATCTGCCTGGG
Reverse Sequence:	GATGGTCCTTCCATCTACCGTG
ACCN:	<u>NM_001130923</u> , <u>NM_001130923.1</u> , <u>NM_001130923.2</u> , <u>BC014879</u> , <u>BC014879.2</u> , <u>BC024281</u> , <u>BC020495</u> , <u>BC040495</u> , <u>BC064551</u> , <u>BC075856</u> , <u>BG177165</u> , <u>BT020037</u> , <u>BX421088</u> , <u>BX432293</u> , <u>NM_001130923.3</u>
UniProt ID:	<u>Q9UNT1</u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)