

Product datasheet for **HP233344**

Archaemetzincin 2 (AMZ2) Human qPCR Primer Pair (NM_016627)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	51321
Forward Sequence:	CGCCCTCTAACCTTTGCCCTA
Reverse Sequence:	TCCTCGTGACTGTGTTCTGGAG
ACCN:	<u>BC050709</u> , <u>NM_016627</u> , <u>NM_016627.1</u> , <u>NM_016627.2</u> , <u>NM_016627.3</u> , <u>NM_016627.4</u> , <u>BC050709.1</u> , <u>BC056271</u> , <u>BC056271.1</u> , <u>BC035935</u> , <u>BE543044</u> , <u>BF978559</u> , <u>BM450564</u> , <u>BQ432895</u>
UniProt ID:	<u>Q86W34</u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)