

Product datasheet for **HP233391**

Acyl CoA Thioesterase 9 (ACOT9) Human qPCR Primer Pair (NM_001037171)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	23597
Forward Sequence:	ATCCA CT CCGCCAAGATGTCTC
Reverse Sequence:	GATGTCTTCCCGACCCAGCTAA
ACCN:	<u>BC012573</u> , <u>NM_001037171</u> , <u>NM_001037171.1</u> , <u>BC136671</u> , <u>BC144360</u> , <u>BC144370</u> , <u>BX648512</u> , <u>NM_001037171.2</u>
UniProt ID:	<u>Q9Y305</u>
Synonyms:	ACATE2; CGI-16; MT-ACT48; MTA CT 48
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)