

Product datasheet for **HP233442**

PHF19 Human qPCR Primer Pair (NM_015651)

Product data:

| | |
|-------------------|--|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 26147 |
| Forward Sequence: | TTCTGCTCCGTGTGTAACCAGG |
| Reverse Sequence: | GGTGGTTGACAAAGGCCAGAATC |
| ACCN: | NM_015651 , NM_015651.1 , NM_015651.2 , BC125076 , BC022374 , BC044224 , BC108663 , BC125077 , BM784783 , BM784790 , BU632411 , BX640713 , NM_015651.3 |
| UniProt ID: | Q5T6S3 |
| Synonyms: | MTF2L1; PCL3; TDRD19B |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions) |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)