

Product datasheet for **HP233445**

Sec8 (EXOC4) Human qPCR Primer Pair (NM_021807)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	60412
Forward Sequence:	CTGCAAACCTGGAGCCAGAAAC
Reverse Sequence:	CGAAGAGGACACTGTTTGGCTG
ACCN:	NM_021807 , NM_021807.1 , NM_021807.2 , NM_021807.3 , BC067263 , BC067263.1 , BC020607 , BC026174 , BE244117 , BE246205 , NM_021807.4
UniProt ID:	Q96A65
Synonyms:	SEC8; SEC8L1; Sec8p
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)