

Product datasheet for **HP233454**

CHID1 Human qPCR Primer Pair (NM_001142676)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	66005
Forward Sequence:	GAGTGTGGTTCTTGAGCATCGC
Reverse Sequence:	AGACCTTGGTGACATCGTAGCC
ACCN:	<u>BC000001</u> , <u>NM_001142676</u> , <u>NM_001142676.1</u> , <u>BC000001.2</u> , <u>BC095409</u> , <u>BC095409.1</u> , <u>BC013642</u> , <u>BI759317</u> , <u>BI818660</u> , <u>BI821281</u> , <u>BN000479</u> , <u>NM_001142676.2</u>
UniProt ID:	<u>Q9BWS9</u>
Synonyms:	GL008; SI-CLP; SICLP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)